

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Caspase-8 Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 11014

**Size:** 0.1 mg affinity purified antibody in phosphate-buffered saline, pH 7.4.

### SPECIFICATION SUMMARY

**Antigen:** Peptide corresponding to aa 266-436 of human caspase-8.

**Host Species:** Chicken

**Stabilizers:** None

**Preservatives:** None. Available upon request.

### SPECIFICITY

This antibody recognizes human caspase-8, a member of the ICE/CED-3 family of cysteine proteases. Caspase-8 appears to be essential for TNF-mediated apoptosis.

### APPLICATIONS

*Immunoblotting:* use at 0.5-2 ug/ml. In immunoblots, a band of approximately 55 kD is detected.

### DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

### STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C. Avoid multiple freeze-thaw cycles.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*