

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Bcl-w Polyclonal Antibody

### ORDERING INFORMATION

Catalog Number	Quantity
11056-25	25 ug
11056-100	100 ug

### SPECIFICATION SUMMARY

**Antigen:** A synthetic peptide based on residues 16-29 of human Bcl-w.

**Host Species:** Rabbit

**Format:** Affinity Purified

**Buffer:** Borate buffered saline, pH 8.2

### SPECIFICITY

This antibody detects a 20 kDa protein on rat brain tissue extract SDS-PAGE immunoblots, corresponding to the molecular weight of human Bcl-w. This antibody reacts with human, monkey, mouse, rat, and bovine Bcl-w.

### APPLICATIONS

*Western Blot:* 10 ug/ml

*Immunohistochemistry:* 20 ug/ml

### STORAGE AND STABILITY

Store frozen product at or below -20°C.

Thawed product may be stored for 2-4 weeks at 4°C. For optimal storage, aliquot and store at -20°C.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*