

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Bromodeoxyuridine Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 12521

Size: 100 ug IgG in 10 mM PBS, pH 7.4, lyophilized.

SPECIFICATION SUMMARY

Antigen: Bromodeoxyuridine-BSA conjugate.

Host Species: Mouse

Clone Number: BU6-4

Subclass: IgG1

Stabilizers: 1% BSA

Preservatives: 0.1% Sodium Azide

STORAGE AND STABILITY

Lyophilized antibody is stable at 4°C for 2 years. Store the reconstituted solution in aliquots at -20°C for one year or at 4°C for 6 months. Avoid repeated freeze-thaw cycles.

SPECIFICITY

This antibody specifically reacts with bromodeoxyuridine.

APPLICATIONS

Immunoassay: 5 ug/ml

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

DILUTION INSTRUCTIONS

Reconstitute in distilled water to the final desired concentration.