

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

BRD2 Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13020

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein between the putative nuclear localization signal and the ET domain (LocusLink ID 6046).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.025 – 0.1 ug/ml.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.