

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

CKII alpha Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13028

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein located within exons 12 and 13 (LocusLink ID 1457).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.3 – 1 ug/ml.

Immunoprecipitation: Use Catalog Number 13030.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.