

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## mTOR Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 13034

**Size:** 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

### SPECIFICATION SUMMARY

**Antigen:** Synthetic peptide representing a portion of the protein encoded in part by exons 5 and 6 (LocusLink ID 1981).

**Host Species:** Rabbit

**Stabilizers:** None

**Preservatives:** 0.1% Sodium Azide

### SPECIFICITY

This antibody reacts with human and rat. Antibody is affinity purified.

### APPLICATIONS

*Western Blot:* Use Catalog Number 13036.

*Immunoprecipitation:* 1-4 ug/mg lysate.

### STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*