

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

STK15 Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13042

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein encoded in exons 5 and 6 (LocusLink ID 6790 and 8465).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.05-0.2 ug/ml.

Immunoprecipitation: 1-5 ug/mg lysate.

Immunofluorescence: 1 ug/ml.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.