

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Epstein Barr Virus Monoclonal Antibodies

ORDERING INFORMATION

Catalog No.: 18501(1H1) and 18502 (4A8)

Size: 100 ug Protein A-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Epstein Barr Virus.

Host Species: Mouse

Antibody Class: 18501 = IgG1;
18502 = IgG2b.

Stabilizers: None

Preservatives: 0.1% sodium azide.

SPECIFICITY

18501 reacts with p120 and p160 capsid proteins and 18502 reacts with p120 capsid protein.

APPLICATIONS

These antibodies may be used in ELISA and immunofluorescence.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

*For in vitro investigational use only.
Not for use in therapeutic or diagnostic procedures.*