

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Hepatitis A Virus Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 18603

Size: 100 ug Protein A-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Purified Hepatitis A virus.

Host Species: Mouse

Mouse Isotype: IgG1

Stabilizers: None

Preservatives: 0.1% sodium azide.

SPECIFICITY

This antibody reacts Hepatitis A virus (epitope unknown).

APPLICATIONS

This antibody may be used in ELISA and immunohistochemistry to detect Hepatitis A virus.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C to -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.