

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

6-His Affinity Purified Antibody, Agarose Immobilized

ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18812A	0.1 mg

Buffer: 0.01M Phosphate, pH 7.2, 0.1M NaCl

Stabilizer: None

Preservatives: 0.1% sodium azide

Concentration: 100 ug antibody / 200 ul gel

SPECIFICATION SUMMARY

Antigen: 6-His (HHHHHH) coupled to KLH. Antibody coupled to agarose beads using a cyanogen bromide method.

Host Species: Goat

APPLICATIONS

Immunoprecipitation: 15 to 25 ul of gell slurry per 0.1 to 1 mg of protein lysate or extract

STORAGE AND STABILITY

Store at 4°C; stable for one year

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.