

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## 6-His Affinity Purified Antibody, HRP Conjugated

### ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18814P	0.1 mg

**Buffer:** 50mM HEPES, pH 7.1, 0.1M NaCl

**Stabilizer:** 0.2% BSA

**Preservative:** 0.1% Pro-Clean 400

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** 6-His (HHHHHH) coupled to KLH.

**Host Species:** Rabbit

### APPLICATIONS

*ELISA:* 0.01 – 0.1 ug/ml

*Western Blot:*

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

*IHC/ICC:* 0.5-2 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*