

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Glutathione-S-Transferase (GST) Affinity Purified Antibody, FITC Conjugated

### ORDERING INFORMATION

<b><u>Catalog No.</u></b>	<b><u>Size:</u></b>
18842F	0.1 mg

**Buffer:** Phosphate Buffered Saline, pH 7.2.

**Stabilizer:** 0.2% BSA

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** GST from *Schistosoma japonicum* coupled to KLH

**Host Species:** Goat

**Preservatives:** 0.1% sodium azide

### APPLICATIONS

*IHC/ICC:* 0.5-5 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*