

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Glutathione-S-Transferase (GST) Affinity Purified Antibody, HRP Conjugated

ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18842P	0.1 mg

Buffer: 50mM HEPES, pH 7.1, 0.1M NaCl

Stabilizer: 0.2% BSA

Preservatives: 0.01% thimerosal

Concentration: 1 mg/ml

SPECIFICATION SUMMARY

Antigen: GST from *Schistosoma japonicum* coupled to KLH.

Host Species: Goat

APPLICATIONS

ELISA: 0.01 – 0.1 ug/ml

Western Blot:

Colorimetric: 0.1-1 ug/ml

Chemiluminescent: 0.03-1 ug/ml

IHC/ICC: 0.5-2 ug/ml

STORAGE AND STABILITY

Store at 4°C; stable for one year

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.