

# QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

## ***Staph. aureus* Enterotoxin Monoclonal Antibodies**

### **ORDERING INFORMATION**

**Catalog No.:** 20401 and 20402

**Size:** 200 ug Protein A-purified antibody.

### **SPECIFICATION SUMMARY**

**Antigen:** Purified enterotoxin A from FRI-722H strain of *Staph. aureus*.

**Host Species:** Mouse

**Antibody Class:** IgG2a

**Stabilizers:** None

**Preservatives:** 0.1% sodium azide

### **DILUTION INSTRUCTIONS**

Dilute in PBS or medium which is identical to that used in the assay system.

### **APPLICATIONS**

These antibodies may be used in immunoassays to detect *Staph. aureus* enterotoxins.

### **STORAGE AND STABILITY**

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

### **REFERENCE**

Advances In Enzyme Regulation (ed. by G. Weber), vol. 30, pp. 331-335, Pergamon Press , New York, 1990.

### **% CROSS-REACTIVITY WITH *Staph. aureus* ENTEROTOXINS**

	<b>A</b>	<b>B</b>	<b>C1</b>	<b>C2</b>	<b>D</b>	<b>E</b>
20401	21	0	0	0	100	16
20402	8	0	0	0	100	0

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*