

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Human Plasminogen Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 21401-21403

Size: 500 ug Protein A-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Purified human plasminogen.

Host Species: Mouse

Antibody Class: 21401 = IgG2b;
21402 = IgG1; 21403 = IgG2a.

Stabilizers: None

Preservatives: 0.1% sodium azide.

SPECIFICITY

These antibodies are specific for human plasminogen.

APPLICATIONS

These antibodies may be used in immunoassays to detect/quantitate human plasminogen. They may also be used to affinity purify plasminogen and in immunohistochemistry on formalin-fixed/paraffin embedded tissues.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.