

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Anti-Apaf1 (NT) Antibody

ORDERING INFORMATION

Catalog No.: 2412

Size: 100 ug IgG in PBS, pH 7.4, purified by immunoaffinity chroma-tography.

BACKGROUND

Apoptosis is related to many diseases and induced by a family of cell death receptors and their ligands. Cell death signals are transduced by death domain containing adapter molecules and members of the caspase family of proteases. The mammalian homologues of the key cell death gene CED-4 in *C. elegans* has been identified recently from human and mouse and designated Apaf1 (for apoptosis protease-activating factor 1). Apaf1 binds to cytochrome c (Apaf2) and caspase-9 (Apaf3), which leads to caspase-9 activation. Activated caspase-9 in turn cleaves and activates caspase-3, one of the key proteases responsible for the proteolytic cleavage of many key proteins in apoptosis. Apaf1 can also associate with caspase-4 and caspase-8. Apaf1 is ubiquitously expressed in human tissues.

SPECIFICATION SUMMARY

Antigen: Peptide corresponding to aa 12-28 of human Apaf1. The sequences of the immunogenic peptide are identical between human and mouse.

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.02% sodium azide.

SPECIFICITY

This antibody recognizes full-length Apaf1 (130 kD), mouse and human.

APPLICATIONS

Immunoblotting: use at 1:1,000-1:2,000 dilution.

Positive control: Whole cell lysate of HeLa cells

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.