

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

MMP 16b Antibody

ORDERING INFORMATION

Catalog No.: 29018

Size: 1 mg purified IgY in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Peptide corresponding to aa 407-458 of human MMP 16b (isoform 2 of MMP 16).

Host Species: Chicken

Preservatives: None. Available upon request.

Stabilizers: None

SPECIFICITY

This antibody specifically recognizes human MMP 16. It does not cross-react with other members of the MMP family.

APPLICATIONS

Immunoblotting: use at 0.5-2 ug/ml. In immunoblots, a band of approximately 70 kD is detected.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C to -70°C. Store antibody in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.