

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

CA-125 Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 3330

Size: 100 ug Protein A-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: CA-125 purified from human ovarian carcinoma.

Host Species: Mouse

Antibody Class: IgG1.

Stabilizers: None

Preservatives: 0.1% sodium azide.

SPECIFICITY

This antibody specifically reacts with human ovarian carcinoma cells

APPLICATIONS

This antibody may be used in immunohistochemistry and flow cytometry to detect CA-125.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C to -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.