

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Human Macrophage Inflammatory-3 Monoclonal Antibody

ORDERING INFORMATION

| Catalog Number | Size |
|----------------|---------|
| 50174-500 | 500 ug |
| 50174-1000 | 1000 ug |

SPECIFICATION SUMMARY

Antigen: recombinant human MIP-3

Host Species: Mouse

Subtype: IgG1

Stabilizers: None

Preservatives: None

APPLICATIONS

ELISA

DILUTION INSTRUCTIONS

Reconstitute with water. Mix gently, washing the sides of the vial. Wait 1 minute before use.

STORAGE AND STABILITY

Store lyophilized product at 4°C. After reconstitution, aliquot and store at -20°C. Stable 2 years lyophilized; one month in solution at 4°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.