

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

HIV-1 gp120 (PND^{mn}) Monoclonal Antibody

ORDERING INFORMATION

Catalog Number	Size
56012-500	500ug
56012-1000	1000ug

Size: 500ug or 1000ug IgG in PBS, pH 7.4, purified by ion exchange chromatography.

STORAGE AND STABILITY

Store lyophilized product at 4°C. After reconstitution, aliquot and store at -20°C.

SPECIFICATION SUMMARY

Antigen: r.gp120 (MN strain)

Host Species: Mouse

Subtype: IgG1

Clone Number: HIV1gp120

Stabilizers: None

Preservatives: None

APPLICATIONS

ELISA

DILUTION INSTRUCTIONS

Reconstitute with water. Mix gently, washing the sides of the vial. Wait 1 minute before use.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

10919 Technology Place, Suite C
Tel. (800) 929-2114 / (858) 675-2405

San Diego, California 92127
Fax (858) 592-1509

www.qedbio.com
info@qedbio.com