

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Lysosome Associated Membrane Protein 2 (LAMP2) Monoclonal Antibody

ORDERING INFORMATION

Catalog Number	Quantity
56272-50	50 ug
56272-200	200 ug

SPECIFICATION SUMMARY

Antigen: Mouse liver lysosomal membrane preparation.

Host Species: Rat

Isotype: IgG

Clone Number: GL2A7

Format: Ascites Fluid

SPECIFICITY

This antibody will label presumptive lysosomes and late endosomes in saponin-permeabilized cells in Immunocytochemistry applications. This antibody reacts with mouse LAMP2.

APPLICATIONS

Immunoprecipitation

Immunocytochemistry

STORAGE AND STABILITY

Store frozen product at or below -20°C.

Thawed product may be stored for 2-4 weeks at 4°C. For optimal storage, aliquot and store at -20°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.