

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Goat anti Mouse IgG, Heavy and Light Chain Specific Antibody

ORDERING INFORMATION

Catalog No: 80000

Size: 2 mg IgG

SPECIFICATION SUMMARY

Antigen: Mouse IgG

Host Species: Goat

Stabilizers: None

Preservatives: None

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, Sterile Liquid

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chains and light chains associated with mouse immunoglobulin.

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from the date of receipt.

For *in vitro* research use only.