

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Goat anti Mouse IgG, Heavy and Light Chain Specific Antibody, Fluorescein Conjugated

ORDERING INFORMATION

Catalog No: 80300

Size: 2.0 mg, lyophilized powder.

SPECIFICATION SUMMARY

Antigen: Mouse IgG

Host Species: Goat

Stabilizers: 15 mg/mL Bovine Serum Albumin

Preservatives: 0.05% Sodium Azide; 0.01% Thimerosal

Buffer: 0.01M Sodium Phosphate, 0.25M NaCL, pH 7.6.

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chains and light chains associated with mouse immunoglobulin.

APPLICATIONS

Most Applications: 1:50 – 1:200

DILUTION INSTRUCTIONS

Reconstitute with 2.0 mL distilled water.

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from the date of reconstitution.

For *in vitro* research use only.