

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Goat anti Human IgA, Alpha Chain Specific Antibody

### ORDERING INFORMATION

**Catalog No:** 87010

**Size:** 2.0 mg, sterile liquid.

### SPECIFICATION SUMMARY

**Antigen:** Human IgA

**Host Species:** Goat

**Buffer:** 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

**Stabilizers:** None

**Preservatives:** None

**Purification:** Immunoaffinity chromatography

### SPECIFICITY

The antibody reacts specifically with the heavy chain of human IgA immunoglobulin.

### STORAGE AND STABILITY

Store product at 2-8°C until opened. Upon opening, this product should be stored at 2-8°C under sterile conditions. Expires 1 year from the date of receipt.

**For *in vitro* research use only.**