

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Goat anti Human IgM, Mu Chain Specific Antibody

ORDERING INFORMATION

Catalog No: 87020

Size: 2.0 mg, sterile liquid.

SPECIFICATION SUMMARY

Antigen: Human IgM

Host Species: Goat

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

Preservatives: None

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chain associated with human immunoglobulin M.

STORAGE AND STABILITY

Store product at 2-8°C until opened. Upon opening, this product should be stored at 2-8°C under sterile conditions. Expires 1 year from the date of receipt.

For *in vitro* research use only.