

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Goat anti Human IgG (H&L), IgA, IgM Antibody

ORDERING INFORMATION

Catalog No: 87030

Size: 2.0 mg, sterile liquid.

SPECIFICATION SUMMARY

Antigen: Human IgG, IgA, and IgM

Host Species: Goat

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

Stabilizers: None

Preservatives: None

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chains of human IgG, IgA, and IgM immunoglobulin, and with light chains associated with most human immunoglobulins.

STORAGE AND STABILITY

Store product at 2-8°C until opened. Upon opening, this product should be stored at 2-8°C under sterile conditions. Expires 1 year from the date of receipt.

For *in vitro* research use only.