

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Goat anti Human IgA, Alpha Chain Specific Antibody, Alkaline Phosphatase Conjugated

### ORDERING INFORMATION

**Catalog No:** 87210

**Size:** 1.0 ml, lyophilized powder.

### SPECIFICATION SUMMARY

**Antigen:** Human IgA

**Host Species:** Goat

**Buffer:** 0.01M Tris-HCl, 0.25M NaCl, pH 8.0

**Stabilizers:** 15 mg/ml Bovine Serum Albumin

**Preservatives:** None

**Purification:** Immunoaffinity chromatography

### SPECIFICITY

The antibody reacts specifically with the heavy chain associated with human immunoglobulin A.

### APPLICATIONS

ELISA and Western Blotting: 1:5,000 – 1:50,000

### DILUTION INSTRUCTIONS

Reconstitute with 1.0 mL distilled water, or a combination of 0.5 ml of distilled water and 0.5 ml of glycerol if stored frozen.

### STORAGE AND STABILITY

Store at 2-8°C until reconstituted. Upon reconstitution, this product should be stored frozen in 50% glycerol. Expires 1 year from the date of reconstitution.

**For *in vitro* research use only.**