

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

Goat anti Human IgM, Mu Chain Specific Antibody, Alkaline Phosphatase Conjugated

ORDERING INFORMATION

Catalog No: 87220

Size: 1.0 ml, lyophilized powder.

SPECIFICATION SUMMARY

Antigen: Human IgM

Host Species: Goat

Buffer: 0.01M Tris-HCl, 0.25M NaCl, pH 8.0

Stabilizers: 15 mg/ml Bovine Serum Albumin

Preservatives: None

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chain associated with human immunoglobulin M.

APPLICATIONS

ELISA and Western Blotting: 1:5,000 – 1:50,000

DILUTION INSTRUCTIONS

Reconstitute with 1.0 mL distilled water, or a combination of 0.5 ml of distilled water and 0.5 ml of glycerol if stored frozen.

STORAGE AND STABILITY

Store at 2-8°C until reconstituted. Upon reconstitution, this product should be stored frozen in 50% glycerol. Expires 1 year from the date of reconstitution.

For *in vitro* research use only.