

# QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

## Goat anti Human IgM, Mu Chain Specific Antibody, FITC Conjugated

### ORDERING INFORMATION

**Catalog No:** 87320

**Size:** 2.0 mg, lyophilized powder.

### SPECIFICATION SUMMARY

**Antigen:** Human IgM

**Host Species:** Goat

**Buffer:** 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

**Stabilizers:** 15 mg/ml Bovine Serum Albumin

**Preservatives:** 0.05% Sodium Azide

**Purification:** Immunoaffinity chromatography

### SPECIFICITY

The antibody reacts with the heavy chains and light chains associated with human immunoglobulin M.

### MOST APPLICATIONS

Begin diluting at: 1:50 – 1:200

### DILUTION INSTRUCTIONS

Reconstitute with 1.5 mL distilled water, or a combination of 0.75 ml of distilled water and 0.75 ml of glycerol if stored frozen.

### STORAGE AND STABILITY

Store at 2-8°C until reconstituted. Upon reconstitution, this product should be stored frozen in 50% glycerol. Expires 1 year from the date of reconstitution.

**For *in vitro* research use only.**