

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

E2F-1 Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 3334 (clone 1G10)

Size: 100 ug Protein G-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: GST fusion protein expressed in E. coli corresponding to aa 162-476 of human E2F-1.

Host Species: Mouse

Antibody Class: IgG1/ κ

Stabilizers: None

Preservatives: None. Available on request.

SPECIFICITY

This antibody recognizes an epitope within aa 244-264 of human E2F-1.

APPLICATIONS

Immunoblotting: use at 1-5 ug/ml.

Immunoprecipitation: use at 1-5 ug/ml.

Positive controls: E2F-1 transfected Soas2 cells and recombinant fusion protein.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.