

# QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

## Hepatitis C Virus Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 18604

**Size:** 100 ug affinity-purified antibody  
in PBS, pH 7.4

### SPECIFICATION SUMMARY

**Antigen:** Peptide corresponding to  
amino acids 9-21 of Hepatitis C Virus  
(HCV) core protein C1.

**Host Species:** Rabbit

**Stabilizers:** None

**Preservatives:** 0.02% sodium azide.

### SPECIFICITY

This antibody recognizes the core  
protein C1 of Hepatitis C Virus.

### APPLICATIONS

*Immunoblotting:* use at 1:500-1:1,000.

### DILUTION INSTRUCTIONS

Dilute in PBS or medium which is  
identical to that used in the assay  
system.

### STORAGE AND STABILITY

This antibody is stable for at least one  
(1) year at -70°C. Avoid multiple freeze-  
thaw cycles.

*For in vitro investigational use only. Not  
for use in therapeutic or diagnostic  
procedures.*