

# QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

## Hepatitis G Virus Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 18607

**Size:** 100 ug affinity-purified antibody  
in PBS, pH 7.4

### SPECIFICATION SUMMARY

**Antigen:** Peptide corresponding to amino acids 2268-2276 of the NS5 non-structural polyprotein of Hepatitis G Virus (HGV).

**Host Species:** Rabbit

**Stabilizers:** None

**Preservatives:** 0.02% sodium azide.

### SPECIFICITY

This antibody recognizes the NS5 polyprotein of HGV, also known as GBV-C/HBV.

### APPLICATIONS

*Immunoblotting:* use at 1:500-1:1,000.

### DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

### STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -70°C. Avoid multiple freeze-thaw cycles.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*