

QED Bioscience Inc.

ADVANCED ANTIBODY TECHNOLOGIES

PRAC Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 56245

Size: 100 ug affinity-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Human PRAC protein.

Host Species: Rabbit

Antibody Class: Polyclonal

Stabilizers: None

Preservatives: 0.02% sodium azide.

SPECIFICITY

This antibody recognizes human PRAC protein, a gene product that is expressed in normal human prostate, rectum, and colon as well as in prostate carcinoma. PRAC is localized in cell nuclei and may play a regulatory role in the nucleus.

APPLICATIONS

Immunoblotting: use at 1:500-1:1,000 dilution. A band of approximately 6 kD is detected.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.