

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Rabbit (IgG) anti Chicken IgY, Heavy and Light Chain Specific Antibody

ORDERING INFORMATION

Catalog No: 86000

Size: 2.0 mg IgG

SPECIFICATION SUMMARY

Antigen: Chicken IgY

Host Species: Rabbit

Stabilizers: None

Preservatives: None

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, Sterile Liquid

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chains and light chains associated with chicken immunoglobulin.

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from the date of receipt.

For *in vitro* research use only.