

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Rabbit (IgG) anti Chicken IgY, Heavy and Light Chain Specific Antibody, Peroxidase Conjugated

ORDERING INFORMATION

Catalog No: 86100

Size: 1.5 mL, lyophilized powder.

SPECIFICATION SUMMARY

Antigen: Chicken IgY

Host Species: Rabbit

Stabilizers: 15 mg/mL Bovine Serum Albumin

Preservatives: 0.01% Thimerosal

Buffer: 0.01M Sodium Phosphate, 0.25M NaCL, pH 7.6.

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chains and light chains associated with chicken immunoglobulin.

APPLICATIONS

Immunohistochemistry: 1:500 – 1:5,000

ELISA or Western Blot, Chromogenic Substrate: 1:5,000 – 1:100,000

Western Blot, ECL Substrate: 1:10,000 – 1:200,000

DILUTION INSTRUCTIONS

Reconstitute with 1.5 mL distilled water.

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from the date of reconstitution.

For *in vitro* research use only.