

Rabbit Polyclonal Antibody Development Services

QED ADVANCED ANTIBODY TECHNOLOGIES

QED Bioscience has over 23 years of experience in the development of antibodies to proteins, peptides, small molecules, and other antigens. Our rabbit polyclonal antibody development program utilizes 2 New Zealand White Rabbits in a 60-day protocol that will provide you with approximately 200 mL of antiserum. ELISA testing of the bleeds against your antigen at the end of the project is included in the base price. All work is done by QED personnel at our facilities.

Your Rabbit Antibody Development Project Includes:

- Two New Zealand White Rabbits
- Project Initiation, Housing and Husbandry, and Termination
- All Immunizations, blood collection, and serum processing
- A 60-day Protocol
- Approximately 180-200 mL of antiserum
- ELISA testing of all bleeds against your antigen at project termination

Base Price for the Above Services: \$675.00

Pricing discounts are available based upon the number of projects initiated over a 12 month period.

Additional Services are Available:

- ELISA Testing at each bleed \$75.00 per bleed
- Affinity Purification \$595.00 for 50 mL: \$2.50 per mL over 50 mL
- Conjugation of Antigen to KLH or BSA \$200.00
- A Complete Peptide Antibody Development Package is Available.
Please contact QED for a customized price quote.

For more information about our products and services please contact us at:

Phone: 800-929-2114 / 858-675-2405

Fax: 858-592-1509

e-mail: info@qedbio.com

or visit our website at www.qedbio.com

