

Recombinant Human Prolactin

ORDERING INFORMATION

Catalog Number	Size
12205P-10	10 ug
12205P-50	50 ug
12205P-1000	1000 ug

Format: Sterile filtered, lyophilized.

Purification: Chromatographic

DESCRIPTION

Recombinant Human Leptin produced in *E. coli* is a single, non-glycosylated, polypeptide chain containing 200 amino acids with a molecular mass of 23,007 daltons. This product has full biological activity.

SPECIFICATION SUMMARY

Source: *Escherichia coli*

Stabilizers: None

Preservatives: None

Purity: Greater than 95% by RP-HPLC, anion-exchange chromatography; reducing and non-reducing SDS-PAGE, silver stain.

Dimers/Aggregates: <1% by SDS-PAGE

Biological Activity: Fully biologically active compared to reference standards. ED₅₀ determined by the dose-dependent proliferation of rat OB-lymphoma Nb2-11 is <0.065 ng/ml.

Endotoxin: <0.1 ng/ug

Protein Content: Determined spectrophotometrically at 280nm and by RP-HPLC analysis against a calibrated standard.

SOLUBILITY

Reconstitute the lyophilized product in sterile water not less than 100 ug/ml. Dialyzed against 10 mM sodium phosphate pH 8 / 50 mM NaCl prior to lyophilization.

STORAGE AND STABILITY

Lyophilized product is stable at room temperature for 3 weeks. Store lyophilized product desiccated at -20°C. Upon reconstitution, store at 4°C for 2-7 days or aliquot and store at -20°C. For long-term storage add either HSA or BSA to 0.1%. Avoid repeated freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures