

c-Jun Polyclonal Antibody

ORDERING INFORMATION

Catalog Number	Size
12514-25	25 ug
12514-100	100 ug

Supplied as affinity purified antibody at 1 mg/ml in borate buffered saline, pH 8.4.

SPECIFICATION SUMMARY

Antigen: Peptide corresponding to aa 266-285 of human c-Jun.

Host Species: Rabbit

Subclass: Polyclonal

Clone Number: Not applicable.

Stabilizers: None

Preservatives: 0.02% Sodium Azide

SPECIFICITY

This antibody recognizes human, monkey, rat, mouse, hamster, rabbit, ovine, porcine, bovine, canine, chicken, guinea pig, *Xenopus*, and *Drosophila c-Jun*.

APPLICATIONS

Western blot: Use at 1 ug/ml. A band of approximately 46 kD is detected on immunoblots and corresponds to the molecular weight of c-Jun.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

Store frozen product at or below -20°C. Thawed product may be stored for 2-4 weeks at 4°C. For optimal storage, aliquot to smaller portions and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.