

## 6-His Affinity Purified Antibody, HRP Conjugated

### ORDERING INFORMATION

<b><u>Catalog No.</u></b>	<b><u>Size:</u></b>
18812P	0.1 ml

**Buffer:** Phosphate buffered saline (PBS) containing 0.2% BSA and 0.1% Pro-Clean 400

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** 6-His (HHHHHH) coupled to KLH. Antibody was isolated by affinity chromatography using the peptide immobilized on solid support then conjugated to horseradish peroxidase (HRP). NOTE: This antibody does not detect all N-terminal 6-His tagged proteins.

**Host Species:** Goat

### APPLICATIONS

*ELISA:* 1:10,000-1:100,000

*Western Blot:*

*Colorimetric:* 1:1,000-1:10,000

*Chemiluminescent:* 1:1,000-1:30,000

*Immunohistochemistry:* not recommended

*Immunocytochemistry:* 1:200-1:500

### STORAGE AND STABILITY

Store at 2-8°C; stable for one year from date of receipt.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*