

## AU5 Affinity Purified Antibody

### ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18820	1 mg
18820-01	100 ug

**Buffer:** Phosphate Buffered Saline, pH 7.2.

**Concentration:** 1 mg/ml

**Preservatives:** 0.1% sodium azide

### SPECIFICATION SUMMARY

**Antigen:** AU5 (TDFYLK) coupled to KLH.

**Host Species:** Goat

### APPLICATIONS

*ELISA:*

*Coating:* 2-10 ug/ml

*Primary:* 0.03-1 ug/ml

*Western Blot:*

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

*IHC/ICC:* 0.5-5 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*