

## c-myc Affinity Purified Antibody, HRP Conjugated

### ORDERING INFORMATION

<b><u>Catalog No.</u></b>	<b><u>Size:</u></b>
18824P	0.1 mg

**Buffer:** 50mM HEPES, pH 7.1, 0.1M NaCl

**Stabilizer:** 0.2% BSA

**Preservatives:** 0.01% thimerosal

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** c-myc (EQKLISEEDL) coupled to KLH.

**Host Species:** Goat

### APPLICATIONS

**ELISA:** 0.01 – 0.1 ug/ml

**Western Blot:**

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

**IHC/ICC:** 0.5-2 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*