

## glu-glu Affinity Purified Antibody, HRP Conjugated

### **ORDERING INFORMATION**

<b><u>Catalog No.</u></b>	<b><u>Size:</u></b>
18836P	0.1 mg

**Buffer:** 50mM HEPES, pH 7.1, 0.1M NaCl

**Stabilizer:** 0.2% BSA

**Preservatives:** 0.01% thimerosal

**Concentration:** 1 mg/ml

### **SPECIFICATION SUMMARY**

**Antigen:** glu-glu (EYMPEME) coupled to KLH.

**Host Species:** Goat

### **APPLICATIONS**

**ELISA:** 0.01 – 0.1 ug/ml

**Western Blot:**

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

**IHC/ICC:** 0.5-2 ug/ml

### **STORAGE AND STABILITY**

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*