

## glu-glu Affinity Purified Antibody, HRP Conjugated

### ORDERING INFORMATION

|                           |                     |
|---------------------------|---------------------|
| <b><u>Catalog No.</u></b> | <b><u>Size:</u></b> |
| 18838P                    | 0.1 mg              |

**Buffer:** 50mM HEPES, pH 7.1, 0.1M NaCl

**Stabilizer:** 0.2% BSA

**Preservatives:** 0.01% thimerosal

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** glu-glu (EYMPME) coupled to KLH.

**Host Species:** Rabbit

### APPLICATIONS

*ELISA:* 0.01 – 0.1 ug/ml

*Western Blot:*

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

*IHC/ICC:* 0.5-2 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*