

## Glutathione-S-Transferase (GST) Affinity Purified Antibody, HRP Conjugated

### ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18844P	0.1 mg

**Buffer:** 50mM HEPES, pH 7.1, 0.1M NaCl

**Stabilizer:** 0.2% BSA

**Preservatives:** 0.1% Pro-Clean 400 added

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** GST from *Schistosoma japonicum* coupled to KLH

**Host Species:** Rabbit

### APPLICATIONS

**ELISA:** 0.01 – 0.1 ug/ml

**Western Blot:**

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

**IHC/ICC:** 0.5-2 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*