

## KT3 Affinity Purified Antibody

### **ORDERING INFORMATION**

<b><u>Catalog No.</u></b>	<b><u>Size:</u></b>
18854	1 mg
18854-01	100 ug

**Buffer:** Phosphate Buffered Saline, pH 7.2.

**Concentration:** 1 mg/ml

**Preservatives:** 0.1% sodium azide

### **SPECIFICATION SUMMARY**

**Antigen:** KT3 (KPPTPPPEPET) coupled to KLH.

**Host Species:** Goat

### **APPLICATIONS**

*ELISA:*

*Coating:* 2-10 ug/ml

*Primary:* 0.03-1 ug/ml

*Western Blot:*

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

*IHC/ICC:* 0.5-5 ug/ml

### **STORAGE AND STABILITY**

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*