

T7 Affinity Purified Antibody

ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18862	1 mg
18862-01	100 ug

Buffer: Phosphate Buffered Saline, pH 7.2.

Concentration: 1 mg/ml

Preservatives: 0.1% sodium azide

SPECIFICATION SUMMARY

Antigen: T7 (MASMTGGQMG) coupled to KLH

Host Species: Goat

APPLICATIONS

ELISA:

Coating: 2-10 ug/ml

Primary: 0.03-1 ug/ml

Western Blot:

Colorimetric: 0.1-1 ug/ml

Chemiluminescent: 0.03-1 ug/ml

IHC/ICC: 0.5-5 ug/ml

STORAGE AND STABILITY

Store at 4°C; stable for one year

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.