

## VSV-G Affinity Purified Antibody

### ORDERING INFORMATION

<u>Catalog No.</u>	<u>Size:</u>
18874	1 mg
18874-01	100 ug

**Buffer:** Phosphate Buffered Saline, pH 7.2.

**Concentration:** 1 mg/ml

### SPECIFICATION SUMMARY

**Antigen:** VSV-G (YTDIEMNRLGK) coupled to KLH

**Host Species:** Goat

**Preservatives:** 0.1% sodium azide

### APPLICATIONS

*ELISA:*

*Coating:* 2-10 ug/ml

*Primary:* 0.03-1 ug/ml

*Western Blot:*

*Colorimetric:* 0.1-1 ug/ml

*Chemiluminescent:* 0.03-1 ug/ml

*IHC/ICC:* 0.5-5 ug/ml

### STORAGE AND STABILITY

Store at 4°C; stable for one year

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*