

Fas Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 2432, 2433

Size: 100 ug Protein A-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Fas antigen (31 kD) purified from Jurkat and HeLa cell lines and conjugated to bovine serum albumin.

Antibody Class: IgG1

Host Species: Mouse

Stabilizers: None

Preservatives: 0.1% sodium azide.

SPECIFICITY

These antibodies recognize two different epitopes on the extracellular fas antigen domain.

APPLICATIONS

This antibody may be used in ELISA and Western blot on cell lysates to detect fas antigen.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C to -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.