

Truncated Retinoblastoma Protein (p56^{RB})

ORDERING INFORMATION

Catalog No.: 3109

Size: 50 ug

Concentration: 5.4 mg/ml.

STORAGE AND STABILITY

This product is stable for at least six (6) months at -70°C or below.

SPECIFICATION SUMMARY

Source: Genetically engineered *E. coli*.
Corresponds to amino acids 379-928.

Molecular Weight: 56,000 daltons.

Purity: >90%.

Diluent: 20 mM sodium phosphate
200 mM sodium chloride
1 mM EDTA
10% glycerol
pH 7.5

For in vitro investigational use only.

Not for use in diagnostic or therapeutic procedures.

Preservatives: None.

Biologically Active: No.

Activity Assay: [³H]-Thymidine incorporation
anti-proliferation bioassay.