

Parathyroid Hormone Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 58004, 58005

Size: 200 ug Protein A-purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Purified human parathyroid hormone.

Host Species: Mouse

Mouse Isotype: IgG1

Stabilizers: None

Preservative: 0.1% sodium azide.

SPECIFICITY

58004 recognizes an epitope located at the C-terminus of human parathyroid hormone (aa 53-83), and 58005 recognizes an epitope located at the N-terminus of human parathyroid hormone (aa 1-34).

APPLICATIONS

These antibodies may be used in ELISA and immunohistochemistry to detect human parathyroid hormone.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C to -70°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.